

# Product Datasheet

## IL-9 Antibody (222622) [PE] FAB4091P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB4091P](http://www.novusbio.com/FAB4091P)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB4091P](http://www.novusbio.com/reviews/destination/FAB4091P)



**FAB4091P**

IL-9 Antibody (222622) [PE]

| Product Information         |                                                                                                                          |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                              |
| Storage                     | Store at 4C in the dark.                                                                                                 |
| Clonality                   | Monoclonal                                                                                                               |
| Clone                       | 222622                                                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                                                       |
| Isotype                     | IgG2b                                                                                                                    |
| Conjugate                   | PE                                                                                                                       |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                               |
| Buffer                      | PBS                                                                                                                      |
| Product Description         |                                                                                                                          |
| Host                        | Rat                                                                                                                      |
| Gene ID                     | 3578                                                                                                                     |
| Gene Symbol                 | IL9                                                                                                                      |
| Species                     | Mouse                                                                                                                    |
| Specificity/Sensitivity     | Detects mouse IL-9 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human IL-9.                  |
| Immunogen                   | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant mouse IL-9; Gln19-Pro144<br>Accession # P15247.1 |
| Product Application Details |                                                                                                                          |
| Applications                | Neutralization                                                                                                           |
| Recommended Dilutions       | Neutralization                                                                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB4091P**

---

|                  |                                |
|------------------|--------------------------------|
| NBP2-34909-100ug | Recombinant Human IL-9 Protein |
| 210-TA-005       | TNF-alpha [Unconjugated]       |
| 409-ML-010       | IL-9 [Unconjugated]            |
| M6000B-1         | IL-6 [HRP]                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB4091P](http://www.novusbio.com/reviews/submit/FAB4091P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

