

# Product Datasheet

## Kallikrein 2 Antibody (426723) [PerCP] FAB4104C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB4104C](http://www.novusbio.com/FAB4104C)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB4104C](http://www.novusbio.com/reviews/destination/FAB4104C)



**FAB4104C****Kallikrein 2 Antibody (426723) [PerCP]**

| <b>Product Information</b>         |                                                                                                                                                                                                                             |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                      |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                 |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                                    |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                                                  |
| <b>Clone</b>                       | 426723                                                                                                                                                                                                                      |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                          |
| <b>Isotype</b>                     | IgG1                                                                                                                                                                                                                        |
| <b>Conjugate</b>                   | PerCP                                                                                                                                                                                                                       |
| <b>Purity</b>                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                  |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                                                         |
| <b>Product Description</b>         |                                                                                                                                                                                                                             |
| <b>Host</b>                        | Mouse                                                                                                                                                                                                                       |
| <b>Gene ID</b>                     | 3817                                                                                                                                                                                                                        |
| <b>Gene Symbol</b>                 | KLK2                                                                                                                                                                                                                        |
| <b>Species</b>                     | Human                                                                                                                                                                                                                       |
| <b>Specificity/Sensitivity</b>     | Detects human Kallikrein 2 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human Kallikrein 1, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, or B1 is observed. |
| <b>Immunogen</b>                   | Mouse myeloma cell line NS0-derived recombinant human Kallikrein 2 Ile25-Pro261<br>Accession # P20151                                                                                                                       |
| <b>Product Application Details</b> |                                                                                                                                                                                                                             |
| <b>Applications</b>                | Western Blot, Immunoprecipitation, Neutralization                                                                                                                                                                           |
| <b>Recommended Dilutions</b>       | Western Blot, Immunoprecipitation, Neutralization                                                                                                                                                                           |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                      |

**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB4104C**

---

|                  |                                            |
|------------------|--------------------------------------------|
| NBP1-97005PCP    | Mouse IgG1 Isotype Control (MG1) [PerCP]   |
| NBP2-51511-0.1mg | Recombinant Human Kallikrein 2 His Protein |
| DKK300           | Kallikrein 3/PSA [HRP]                     |
| 4104-SE-010      | Kallikrein 2 [Unconjugated]                |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB4104C](http://www.novusbio.com/reviews/submit/FAB4104C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

