

# Product Datasheet

## **TNF RI/TNFRSF1A Antibody (55R170) [Dylight 650] FAB430W**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB430W](http://www.novusbio.com/FAB430W)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB430W](http://www.novusbio.com/reviews/destination/FAB430W)



**FAB430W**

TNF RI/TNFRSF1A Antibody (55R170) [Dylight 650]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | 55R170                                                                                                                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                             |
| Conjugate                   | DyLight 650                                                                                                                                                                                                                                                                                                                                                                                     |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                                                      |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Hamster                                                                                                                                                                                                                                                                                                                                                                                         |
| Gene ID                     | 7132                                                                                                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | TNFRSF1A                                                                                                                                                                                                                                                                                                                                                                                        |
| Species                     | Mouse                                                                                                                                                                                                                                                                                                                                                                                           |
| Specificity/Sensitivity     | Detects mouse TNF RI/TNFRSF1A in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human (rh) 4-1BB, rhCD27, rhCD30, recombinant mouse (rm) CD27, rmCD30, rhCD40, rmCD40, rhDR3, rhDR6, rmEDAR, rhFas, rmFas, rhEDAR, rhGITR, rmGITR, rhHVEM, rhLTβR, rmLTβR, rhNGF R, rhOPG, rmOPG, rhRANK, rmRANK, rhTAJ, rhTNF RI, rhTNF RII or rmTNF RII is observed. |
| Immunogen                   | Recombinant mouse TNF RI/TNFRSF1A Extracellular domain                                                                                                                                                                                                                                                                                                                                          |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                               |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Immunoprecipitation, Neutralization                                                                                                                                                                                                                                                                                                                                               |
| Recommended Dilutions       | Western Blot, Immunoprecipitation, Neutralization                                                                                                                                                                                                                                                                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB430W](http://www.novusbio.com/reviews/submit/FAB430W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

