

Product Datasheet

Fc gamma RIII (CD16) Antibody (1001049) [Allophycocyanin] FAB4325A

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB4325A

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB4325A



FAB4325A

Fc gamma RIII (CD16) Antibody (1001049) [Allophycocyanin]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	1001049
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	Allophycocyanin
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	2214
Species	Human
Specificity/Sensitivity	Detects human Fc gamma RIIIA/CD16a in direct ELISAs. In direct ELISAs, 60% cross-reactivity with Fc gamma RIIIB/CD16b is observed.
Immunogen	Human embryonic kidney cell, HEK293-derived human Fc gamma RIIIA/CD16a Gly17-Gln208 Accession # P08637

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready, Immunocytochemistry
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

Fc gamma RIII (CD16) Antibody (1001049) [Allophycocyanin] [FAB4325A] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.



Allophycocyanin (APC)

LASER (nm)	FILTER
Red (633,640)	660/10
EXCITATION MAX (nm)	EMISSION MAX (nm)
650	660



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB4325A

IC003A	Mouse IgG2a Isotype Control (20102) [Allophycocyanin]
210-TA-005	TNF-alpha [Unconjugated]
1960-FC-050	Fc gamma RIII (CD16) [Unconjugated]
202-IL-010	IL-2 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB4325A

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

