

Product Datasheet

Jak3 Antibody (452524) [DyLight 650] FAB46991W

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB46991W

Updated 8/1/2024 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB46991W



FAB46991W

Jak3 Antibody (452524) [DyLight 650]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	452524
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	DyLight 650
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	3718
Gene Symbol	JAK3
Species	Human
Specificity/Sensitivity	Detects human Jak3 in direct ELISAs.
Immunogen	<i>E. coli</i> -derived recombinant human Jak3 Gly46-Thr209 Accession # P52333
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	CyTOF-ready, Immunocytochemistry, Intracellular Staining by Flow Cytometry
Recommended Dilutions	Intracellular Staining by Flow Cytometry, Immunocytochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB46991W

NBP2-27231C	Mouse IgG2b Isotype Control (MPC-11) [DyLight 650]
H00003718-P01-2ug	Recombinant Human Jak3 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBL1-12094	Jak3 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB46991W

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

