

Product Datasheet

S100A2 Antibody (674035) [DyLight 650] FAB4870W

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB4870W

Updated 5/8/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB4870W



FAB4870W

S100A2 Antibody (674035) [DyLight 650]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	674035
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	DyLight 650
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	6273
Gene Symbol	S100A2
Species	Human
Specificity/Sensitivity	Detects human S100A2 in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human (rh) S100A1, A4, A6, A7, A8, A9, A10, A11, A13, A16, rhS100B, rhS100P, recombinant mouse S100A1, A4, A6, A8, A9, A10, A13, or A16 is observed.
Immunogen	<i>E. coli</i> -derived recombinant human S100A2 Met1-Pro97 Accession # P29034
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB4870W

NBP2-27231C	Mouse IgG2b Isotype Control (MPC-11) [DyLight 650]
NBP1-51023-0.1mg	Recombinant Human S100A2 His Protein
H00006273-T01	S100A2 293T Cell Transient Overexpression Lysate
NB100-56559	S100A7/Psoriasin Antibody (47C1068) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB4870W

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

