

# Product Datasheet

## TAFA4/FAM19A4 Antibody (480103) [Allophycocyanin] FAB5099A

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB5099A](http://www.novusbio.com/FAB5099A)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB5099A](http://www.novusbio.com/reviews/destination/FAB5099A)



**FAB5099A**

TAF4/FAM19A4 Antibody (480103) [Allophycocyanin]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 480103  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG1  |
| Conjugate           | Allophycocyanin   |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS   |

| Product Description     |   |
|-------------------------|---|
| Host                    | Mouse   |
| Gene ID                 | 151647  |
| Gene Symbol             | FAM19A4   |
| Species                 | Human   |
| Specificity/Sensitivity | Detects human TAF4/FAM19A4 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human TAF1, 2, 3, or 5 is observed. |
| Immunogen               | <i>E. coli</i> -derived recombinant human TAF4<br>Ser35-Arg140<br>Accession # Q96LR4  |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, CyTOF-ready, Intracellular Staining by Flow Cytometry    |
| Recommended Dilutions       | Western Blot, Intracellular Staining by Flow Cytometry, CyTOF-ready    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB5099A**

---

|               |  |
|---------------|--|
| NBP1-97005APC | Mouse IgG1 Isotype Control (MG1) [Allorhynchocyanin] |
| 5099-TA-050   | TAF4/FAM19A4 [Unconjugated]                          |
| NB200-541     | Complement C4b/d Antibody (16D2)                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB5099A](http://www.novusbio.com/reviews/submit/FAB5099A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

