

# Product Datasheet

## Carboxypeptidase B2/CPB2 Antibody (650801) [HRP] FAB6036H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6036H](http://www.novusbio.com/FAB6036H)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6036H](http://www.novusbio.com/reviews/destination/FAB6036H)



**FAB6036H**

Carboxypeptidase B2/CPB2 Antibody (650801) [HRP]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                               |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Clone                       | 650801                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Preservative                | No Preservative                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Conjugate                   | HRP                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                                                                                             |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Gene ID                     | 1361                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Gene Symbol                 | CPB2                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Specificity/Sensitivity     | Detects human Carboxypeptidase B2/CPB2 in direct ELISAs and Western blots. In Western blots, approximately 20% cross-reactivity with recombinant human (rh) Carboxypeptidase A1, A4, recombinant mouse (rm) Carboxypeptidase A1, A4, B1, and 5% cross-reactivity with rhCarboxypeptidase B1, B2, and E is observed under reducing conditions. No cross-reactivity with any of these proteins is observed under nonreducing conditions. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Carboxypeptidase B2/CPB2<br>Phe23-Val423 (predicted)<br>Accession # Q96IY4                                                                                                                                                                                                                                                                                                       |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                 |





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to FAB6036H

|                    |                                                                 |
|--------------------|-----------------------------------------------------------------|
| NBP1-97074         | Mouse IgG1 Isotype Control (MG1) [HRP]                          |
| H00001361-P01-10ug | Recombinant Human Carboxypeptidase B2/CPB2 GST (N-Term) Protein |
| DY1707             | C-Reactive Protein/CRP [Biotin]                                 |
| NBP3-06961         | Human Carboxypeptidase B2/CPB2 ELISA Kit (Colorimetric)         |

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6036H](http://www.novusbio.com/reviews/submit/FAB6036H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

