

# Product Datasheet

## **BMP-10 Antibody (654319) [Allophycocyanin] FAB6038A**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6038A](http://www.novusbio.com/FAB6038A)

Updated 5/10/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6038A](http://www.novusbio.com/reviews/destination/FAB6038A)



**FAB6038A**

BMP-10 Antibody (654319) [Allophycocyanin]

| Product Information         |                                                                                                                                                                                                                                                                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                 |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                    |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                  |
| Clone                       | 654319                                                                                                                                                                                                                                                                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                          |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                                                                       |
| Conjugate                   | Allophycocyanin                                                                                                                                                                                                                                                                                                                             |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                  |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                                                                                                             |
| Host                        | Rat                                                                                                                                                                                                                                                                                                                                         |
| Gene ID                     | 27302                                                                                                                                                                                                                                                                                                                                       |
| Gene Symbol                 | BMP10                                                                                                                                                                                                                                                                                                                                       |
| Species                     | Mouse                                                                                                                                                                                                                                                                                                                                       |
| Specificity/Sensitivity     | Detects mouse BMP-10 in direct ELISAs. In direct ELISAs, 100% cross-reactivity with recombinant human (rh) BMP-10 and less than 10% cross-reactivity with rhBMP-8a and recombinant mouse (rm) BMP-8b is observed. No cross-reactivity with rmBMP-3, -4, -5, -6, -7, -9, -15, rhBMP-1.1, -2, -3b, or the propeptide of rhBMP-10 is observed. |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human BMP-10 Asn314-Arg421<br>Accession # Q9R229                                                                                                                                                                                                                                    |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                             |
| Applications                | Immunohistochemistry                                                                                                                                                                                                                                                                                                                        |
| Recommended Dilutions       | Immunohistochemistry                                                                                                                                                                                                                                                                                                                        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6038A**

---

|             |                                                     |
|-------------|-----------------------------------------------------|
| IC006A      | Rat IgG2a Isotype Control (54447) [Allophycocyanin] |
| 314-BP-010  | BMP-4 [Unconjugated]                                |
| 2926-BP-025 | BMP-10 [Unconjugated]                               |
| 355-BM-010  | BMP-2 [Unconjugated]                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6038A](http://www.novusbio.com/reviews/submit/FAB6038A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

