

# Product Datasheet

## ROR gamma t/RORC2/NR1F3 Antibody (1181A) [FITC] FAB61093F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB61093F](http://www.novusbio.com/FAB61093F)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB61093F](http://www.novusbio.com/reviews/destination/FAB61093F)



**FAB61093F**

ROR gamma t/RORC2/NR1F3 Antibody (1181A) [FITC]

**Product Information**

<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	1181A
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Conjugate</b>	FITC
<b>Purity</b>	Protein A or G purified from cell culture supernatant
<b>Buffer</b>	PBS

**Product Description**

<b>Host</b>	Rabbit
<b>Gene ID</b>	6097
<b>Species</b>	Human
<b>Specificity/Sensitivity</b>	Detects recombinant human ROR $\gamma$ t/RORC2/NR1F3 by direct ELISAs.
<b>Immunogen</b>	<i>E.coli</i> derived recombinant human ROR $\gamma$ T/RORC2/NR1F3 Met1-Gln100 Accession # P51449

**Product Application Details**

<b>Applications</b>	Intracellular Staining by Flow Cytometry
<b>Recommended Dilutions</b>	Intracellular Staining by Flow Cytometry
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.

**Images**

ROR gamma t/RORC2/NR1F3 Antibody (1181A) [FITC] [FAB61093F] - Vial of FITC conjugated antibody. FITC is optimally excited at 498 nm by the Blue laser (488 nm) and has an emission maximum of 519 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB61093F**

---

NBP2-24892	Rabbit IgG Isotype Control [FITC]
M6000B-1	IL-6 [HRP]
DY421-05	IL-17/IL-17A [Biotin]
NB100-39002	FoxP3 Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB61093F](http://www.novusbio.com/reviews/submit/FAB61093F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

