

# Product Datasheet

## TRANCE/TNFSF11/RANK L Antibody (685809) [HRP] FAB6263H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6263H](http://www.novusbio.com/FAB6263H)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6263H](http://www.novusbio.com/reviews/destination/FAB6263H)



**FAB6263H**

TRANCE/TNFSF11/RANK L Antibody (685809) [HRP]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                  |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Clone                       | 685809                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Preservative                | No Preservative                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG3                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Conjugate                   | HRP                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                                                                                                   |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Gene ID                     | 8600                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Gene Symbol                 | TNFSF11                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Species                     | Human, Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Specificity/Sensitivity     | Detects human TRANCE/TNFSF11/RANK L in direct ELISAs and Western blots. In direct ELISAs and Western blots, 100% cross-reactivity with recombinant mouse (rm) TRANCE and approximately 25% cross-reactivity with recombinant human (rh) LIGHT is observed. No cross-reactivity with rh4-1BB Ligand, rhAPRIL, rhBAFF, rhEDA, rhEDA-A2, rhFas Ligand, rhGITR Ligand, rhOX40 Ligand, rhTNF- $\alpha$ , rhTRAIL, rhTWEAK, or rhVEGI is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human TRANCE/TNFSF11/RANK L Gly64-Asp245<br>Accession # AAC51762                                                                                                                                                                                                                                                                                                                             |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                       |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6263H**

---

|                |                                                 |
|----------------|-------------------------------------------------|
| NBP1-97077     | Mouse IgG3 Isotype Control (MG3) [HRP]          |
| NB100-56392PEP | TRANCE/TNFSF11/RANK L Antibody Blocking Peptide |
| 210-TA-005     | TNF-alpha [Unconjugated]                        |
| 462-TEC-010    | TRANCE/TNFSF11/RANK L [Unconjugated]            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6263H](http://www.novusbio.com/reviews/submit/FAB6263H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

