

Product Datasheet

DNMT3A Antibody (681615) [PerCP] FAB63151C

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB63151C

Updated 11/11/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB63151C



FAB63151C

DNMT3A Antibody (681615) [PerCP]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	681615
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	PerCP
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	1788
Gene Symbol	DNMT3A
Species	Human, Mouse
Specificity/Sensitivity	Detects human DNMT3A in direct ELISAs. In direct ELISAs and Western blots, no cross-reactivity with recombinant human DNMT1 is observed.
Immunogen	<i>E. coli</i> -derived recombinant human DNMT3A Ala353-Lys486 Accession # Q9Y6K1

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB63151C

NBP2-27231PCP	Mouse IgG2b Isotype Control (MPC-11) [PerCP]
H00001788-Q01-10ug	Recombinant Human DNMT3A GST (N-Term) Protein
NBL1-09970	DNMT3A Overexpression Lysate
NB100-56519	DNMT1 Antibody (60B1220.1) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB63151C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

