

Product Datasheet

DNMT3A Antibody (681630) [DyLight 680] FAB6315Y

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB6315Y

Updated 11/11/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB6315Y



FAB6315Y

DNMT3A Antibody (681630) [DyLight 680]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	681630
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 680
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	1788
Gene Symbol	DNMT3A
Species	Human
Specificity/Sensitivity	Detects human DNMT3A in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human DNMT1 is observed.
Immunogen	<i>E. coli</i> -derived recombinant human DNMT3A Ala353-Lys486 Accession # Q9Y6K1
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Immunohistochemistry, Immunocytochemistry
Recommended Dilutions	Immunohistochemistry, Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB6315Y

NBP1-97005FR	Mouse IgG1 Isotype Control (MG1) [DyLight 680]
H00001788-Q01-10ug	Recombinant Human DNMT3A GST (N-Term) Protein
NBL1-09970	DNMT3A Overexpression Lysate
NB100-56519	DNMT1 Antibody (60B1220.1) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB6315Y

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

