

# Product Datasheet

## TLR10 Antibody (670719) [Biotin] FAB6619B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6619B](http://www.novusbio.com/FAB6619B)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6619B](http://www.novusbio.com/reviews/destination/FAB6619B)



**FAB6619B**

TLR10 Antibody (670719) [Biotin]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                       |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                          |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                        |
| Clone                       | 670719                                                                                                                                                                                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                                                                                                             |
| Conjugate                   | Biotin                                                                                                                                                                                                                                                                                                                                            |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                        |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                               |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                   |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                             |
| Gene ID                     | 81793                                                                                                                                                                                                                                                                                                                                             |
| Gene Symbol                 | TLR10                                                                                                                                                                                                                                                                                                                                             |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                             |
| Specificity/Sensitivity     | Detects human TLR10 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) TLR1, 2, 3, 4, 5, 7, 8, 9, recombinant mouse (rm) TLR6, 11, 12, rhIL-1 RI, rhIL-1 RII, rhIL-1 RAcP, rhST2, rhIL-18 R, rhIL-1 Rrp2, rhIL-18 R $\beta$ , rhSIGIRR, rhIL-1 RAPL1, rhTIGIRR, rhMD-1, rhMD-2, rhTIRAP, or rmRP105 is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant human TLR10<br>Asn168-Lys383<br>Accession # Q9BXR5                                                                                                                                                                                                                                                            |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                   |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6619B**

---

|               |                                               |
|---------------|-----------------------------------------------|
| NBP2-27231B   | Mouse IgG2b Isotype Control (MPC-11) [Biotin] |
| NBP3-17185PEP | TLR10 Recombinant Protein Antigen             |
| 210-TA-005    | TNF-alpha [Unconjugated]                      |
| 6619-TR-050   | TLR10                                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6619B](http://www.novusbio.com/reviews/submit/FAB6619B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

