

# Product Datasheet

## **TEX19.1 Antibody (678803) [Allophycocyanin] FAB6645A**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6645A](http://www.novusbio.com/FAB6645A)

Updated 5/10/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6645A](http://www.novusbio.com/reviews/destination/FAB6645A)



**FAB6645A**

TEX19.1 Antibody (678803) [Allophycocyanin]

| Product Information         |                                                                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                           |
| Storage                     | Store at 4C in the dark.                                                                                                                                                              |
| Clonality                   | Monoclonal                                                                                                                                                                            |
| Clone                       | 678803                                                                                                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                    |
| Isotype                     | IgG2a                                                                                                                                                                                 |
| Conjugate                   | Allophycocyanin                                                                                                                                                                       |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                            |
| Buffer                      | PBS                                                                                                                                                                                   |
| Product Description         |                                                                                                                                                                                       |
| Host                        | Rat                                                                                                                                                                                   |
| Gene ID                     | 400629                                                                                                                                                                                |
| Gene Symbol                 | Tex19.1                                                                                                                                                                               |
| Species                     | Mouse                                                                                                                                                                                 |
| Specificity/Sensitivity     | Detects mouse TEX19.1 in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human TEX19, recombinant mouse (rm) TEX19.2, or rmTEX11 is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant mouse TEX19.1<br>Arg189-Pro351<br>Accession # Q99MV2                                                                                              |
| Product Application Details |                                                                                                                                                                                       |
| Applications                | Immunocytochemistry                                                                                                                                                                   |
| Recommended Dilutions       | Immunocytochemistry                                                                                                                                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB6645A**

---

|            |                                                     |
|------------|-----------------------------------------------------|
| IC006A     | Rat IgG2a Isotype Control (54447) [Allophycocyanin] |
| AF1759     | Oct-3/4 Antibody [Unconjugated]                     |
| NBP1-45243 | UBR2 Antibody                                       |
| NB100-2437 | DAZL Antibody                                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6645A](http://www.novusbio.com/reviews/submit/FAB6645A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

