

# Product Datasheet

## **RANK/TNFRSF11A Antibody (80710) [PerCP] FAB6832C**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6832C](http://www.novusbio.com/FAB6832C)

Updated 11/4/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6832C](http://www.novusbio.com/reviews/destination/FAB6832C)



**FAB6832C****RANK/TNFRSF11A Antibody (80710) [PerCP]**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                             |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                                                                                                                      |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                 |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                    |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                                                                                                                                                  |
| <b>Clone</b>                       | 80710                                                                                                                                                                                                                                                                                                                       |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                          |
| <b>Isotype</b>                     | IgG2a                                                                                                                                                                                                                                                                                                                       |
| <b>Conjugate</b>                   | PerCP                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                  |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                                                                                                                                                         |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                             |
| <b>Host</b>                        | Mouse                                                                                                                                                                                                                                                                                                                       |
| <b>Gene ID</b>                     | 8792                                                                                                                                                                                                                                                                                                                        |
| <b>Species</b>                     | Human                                                                                                                                                                                                                                                                                                                       |
| <b>Specificity/Sensitivity</b>     | Detects human RANK/TNFRSF11A in ELISA. In direct ELISAs, no cross-reactivity with recombinant human (rh) 4-1BB, rhBAFF R, rhCD27, rhCD30, rhCD40, rhDR3, rhDR6, rhEDAR, rhFas, rhGITR, rhHVEM, rhLymphotoxin $\beta$ R, rhNGF R, rhOPG, rhTAJ, rhTNF RI, rhTNF RII, or recombinant mouse (rm) OX40, and rmRANK is observed. |
| <b>Immunogen</b>                   | Mouse myeloma cell line NS0-derived recombinant human RANK/TNFRSF11A extracellular domain<br>Accession # Q9Y6Q6                                                                                                                                                                                                             |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                             |
| <b>Applications</b>                | Immunohistochemistry, Immunocytochemistry                                                                                                                                                                                                                                                                                   |
| <b>Recommended Dilutions</b>       | Immunohistochemistry, Immunocytochemistry                                                                                                                                                                                                                                                                                   |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                      |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6832C**

---

|                |                                          |
|----------------|------------------------------------------|
| NB100-56396PEP | RANK/TNFRSF11A Antibody Blocking Peptide |
| 210-TA-005     | TNF-alpha [Unconjugated]                 |
| 683-RK-100     | RANK/TNFRSF11A [Unconjugated]            |
| M6000B-1       | IL-6 [HRP]                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6832C](http://www.novusbio.com/reviews/submit/FAB6832C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

