

# Product Datasheet

## CD23/Fc epsilon RII Antibody (691632) [DyLight 488] FAB6900K

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6900K](http://www.novusbio.com/FAB6900K)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6900K](http://www.novusbio.com/reviews/destination/FAB6900K)



**FAB6900K**

CD23/Fc epsilon RII Antibody (691632) [DyLight 488]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 691632  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG2a   |
| Conjugate           | DyLight 488   |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate  |

| Product Description     |  |
|-------------------------|--|
| Host                    | Rat  |
| Gene ID                 | 2208   |
| Gene Symbol             | Fcer2a   |
| Species                 | Mouse  |
| Specificity/Sensitivity | Detects mouse CD23/Fcε RII in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) Fcε RIA or rhFcε RIG is observed. |
| Immunogen               | Mouse myeloma cell line NS0-derived recombinant mouse CD23/Fcε RII<br>Glu50-Pro331<br>Accession # P20693   |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Flow Cytometry, CyTOF-ready  |
| Recommended Dilutions       | Flow Cytometry, CyTOF-ready  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB6900K**

---

|                |   |
|----------------|---|
| NBP2-35242-5ug | Recombinant Human CD23/Fc epsilon RII Protein |
| 210-TA-005     | TNF-alpha [Unconjugated]                      |
| DCD230         | CD23/Fc epsilon RII [HRP]                     |
| M6000B-1       | IL-6 [HRP]                                    |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6900K](http://www.novusbio.com/reviews/submit/FAB6900K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

