

Product Datasheet

CD155/PVR Antibody (690912) [DyLight 594] FAB6909M

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB6909M

Updated 8/1/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB6909M



FAB6909M

CD155/PVR Antibody (690912) [DyLight 594]

| Product Information | |
|----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 690912 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | DyLight 594 |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|--------------------------------|---|
| Host | Rat |
| Gene ID | 5817 |
| Gene Symbol | PVR |
| Species | Mouse |
| Specificity/Sensitivity | Detects mouse CD155/PVR in direct ELISAs. |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant mouse CD155/PVR Asp29-Leu348 Accession # NP_081790 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|------------------------------------|---|
| Applications | Western Blot, Flow Cytometry, Blockade of Receptor-ligand Interaction, CyTOF-ready, Immunocytochemistry |
| Recommended Dilutions | Western Blot, Flow Cytometry, Blockade of Receptor-ligand Interaction, Immunocytochemistry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB6909M

| | |
|---------------|---|
| NBP3-21304PEP | CD155/PVR Recombinant Protein Antigen |
| DKK300 | Kallikrein 3/PSA [HRP] |
| 9174-CD-050 | CD155/PVR [Unconjugated] |
| MAB666-100 | DNAM-1/CD226 Antibody (102511) [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB6909M

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

