

Product Datasheet

CXCR6 Antibody (1010910) [Allophycocyanin] FAB6991A

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB6991A

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB6991A



FAB6991A**CXCR6 Antibody (1010910) [Allophycocyanin]****Product Information**

| | |
|----------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 1010910 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | Allophycocyanin |
| Purity | Protein A or G purified from ascites |
| Buffer | PBS |

Product Description

| | |
|--------------------------------|--|
| Host | Mouse |
| Gene ID | 10663 |
| Species | Human |
| Specificity/Sensitivity | Detects human CXCR6 in direct ELISAs. |
| Immunogen | HEK293 human embryonic kidney cell line transfected with human CXCR6 Met1-Leu342 Accession # O00574 |

Product Application Details

| | |
|------------------------------|--|
| Applications | Flow Cytometry, CyTOF-ready |
| Recommended Dilutions | Flow Cytometry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

CXCR6 Antibody (1010910) [Allophycocyanin] [FAB6991A] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB6991A

| | |
|---------------|---|
| IC003A | Mouse IgG2a Isotype Control (20102) [Allophycocyanin] |
| NBP1-76871PEP | CXCR6 Antibody Blocking Peptide |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| NBP3-27473 | Human CXCR6 ELISA Kit (Colorimetric) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB6991A

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

