

# Product Datasheet

## Desmocollin-3 Antibody (765129) [PE] FAB7265P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB7265P](http://www.novusbio.com/FAB7265P)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB7265P](http://www.novusbio.com/reviews/destination/FAB7265P)



**FAB7265P**

Desmocollin-3 Antibody (765129) [PE]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 765129                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2a                                                                                       |
| Conjugate           | PE                                                                                          |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                       |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rat                                                                                                                                                                                                                                                                                   |
| Gene ID                 | 1825                                                                                                                                                                                                                                                                                  |
| Gene Symbol             | DSC3                                                                                                                                                                                                                                                                                  |
| Species                 | Mouse                                                                                                                                                                                                                                                                                 |
| Specificity/Sensitivity | Detects mouse Desmocollin-3 in direct ELISAs and Western blots. In direct ELISAs and Western blots, less than 20% cross-reactivity with recombinant human Desmocollin-3 is observed and no cross-reactivity with recombinant mouse (rm) Desmocollin-1 or rmDesmocollin-2 is observed. |
| Immunogen               | Mouse myeloma cell line NS0-derived recombinant mouse Desmocollin-3 Arg136-Leu695<br>Accession # P55850                                                                                                                                                                               |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry                                     |
| Recommended Dilutions       | Western Blot, Immunohistochemistry                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB7265P**

---

|               |                                           |
|---------------|-------------------------------------------|
| NBP2-56621PEP | Desmocollin-3 Recombinant Protein Antigen |
| 664-LI-025    | LIGHT/TNFSF14 [Unconjugated]              |
| 4630-DC-050   | Desmocollin-3                             |
| AF743         | VEGFR3/Flt-4 Antibody [Unconjugated]      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB7265P](http://www.novusbio.com/reviews/submit/FAB7265P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

