

Product Datasheet

TL1A/TNFSF15 Antibody (2116A) [DyLight 550] FAB74422L

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB74422L

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB74422L



FAB74422L

TL1A/TNFSF15 Antibody (2116A) [DyLight 550]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	2116A
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 550
Purity	Protein A or G purified from cell culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	9966
Species	Human
Specificity/Sensitivity	Detects human TL1A/TNFSF15 in direct ELISAs.
Immunogen	<i>E. coli</i> -derived recombinant human TL1A/TNFSF15 Leu72-Leu251 Accession # O95150
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready, Neutralization
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, Neutralization, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB74422L

NBP2-24891R	Rabbit IgG Isotype Control [DyLight 550]
NBP2-35067-5ug	Recombinant Human TL1A/TNFSF15 Protein
210-TA-005	TNF-alpha [Unconjugated]
1896-TL-010	TL1A/TNFSF15 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB74422L

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

