

Product Datasheet

Integrin alpha 2b beta 3 Antibody (2530A) [Allophycocyanin/Cy7] **FAB76161APCCY7**

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB76161APCCY7

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB76161APCCY7



FAB76161APCCY7

Integrin alpha 2b beta 3 Antibody (2530A) [Allophycocyanin/Cy7]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Monoclonal
Clone	2530A
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	Allophycocyanin/Cy7
Purity	Protein A or G purified from cell culture supernatant
Buffer	PBS
Product Description	
Host	Rabbit
Gene ID	3674
Species	Human
Specificity/Sensitivity	Detects Human Integrin $\alpha 2b\beta 3$ heterodimer in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) Integrin $\alpha 2b$, rhIntegrin $\alpha 5$, rhIntegrin $\alpha 5\beta 6$, rhIntegrin $\alpha 8\beta 1$, rhIntegrin $\beta 1$, rhIntegrin $\beta 2$, rhIntegrin $\beta 3$, rhIntegrin $\beta 5$, rhIntegrin $\beta 7$, and recombinant mouse Integrin $\alpha 2b\beta 3$.
Immunogen	Chinese Hamster Ovary cell line, CHO-derived heterodimer of Human Integrin $\alpha 2b$ (Leu32-Arg993, Accession P08514) and Human Integrin $\beta 3$ (Gly27-Asp718, Accession P05106)
Product Application Details	
Applications	Flow Cytometry
Recommended Dilutions	Flow Cytometry
Application Notes	Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB76161APCCY7

NBP2-24891APCCY7	Rabbit IgG Isotype Control [Allophycocyanin/Cy7]
210-TA-005	TNF-alpha [Unconjugated]
7148-A2-025	Integrin alpha 2b beta 3 [Unconjugated]
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB76161APCCY7

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

