

Product Datasheet

Integrin alpha 2b beta 3 Antibody (2530A) [DyLight 550] FAB76161L

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB76161L

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB76161L



FAB76161L

Integrin alpha 2b beta 3 Antibody (2530A) [DyLight 550]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	2530A
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 550
Purity	Protein A or G purified from cell culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	3674
Species	Human
Specificity/Sensitivity	Detects Human Integrin $\alpha 2\beta 3$ heterodimer in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) Integrin $\alpha 2\beta$, rhIntegrin $\alpha 5$, rhIntegrin $\alpha 5\beta 6$, rhIntegrin $\alpha 8\beta 1$, rhIntegrin $\beta 1$, rhIntegrin $\beta 2$, rhIntegrin $\beta 3$, rhIntegrin $\beta 5$, rhIntegrin $\beta 7$, and recombinant mouse Integrin $\alpha 2\beta 3$.
Immunogen	Chinese Hamster Ovary cell line, CHO-derived heterodimer of Human Integrin $\alpha 2\beta$ (Leu32-Arg993, Accession P08514) and Human Integrin $\beta 3$ (Gly27-Asp718, Accession P05106)
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Flow Cytometry, CyTOF-ready
Recommended Dilutions	Flow Cytometry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB76161L

NBP2-24891R	Rabbit IgG Isotype Control [DyLight 550]
210-TA-005	TNF-alpha [Unconjugated]
7148-A2-025	Integrin alpha 2b beta 3 [Unconjugated]
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB76161L

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

