

# Product Datasheet

## **Histone Deacetylase 2/HDAC2 Antibody (806116) [DyLight 405] FAB7679E**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB7679E](http://www.novusbio.com/FAB7679E)

Updated 5/10/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB7679E](http://www.novusbio.com/reviews/destination/FAB7679E)



**FAB7679E**

Histone Deacetylase 2/HDAC2 Antibody (806116) [DyLight 405]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 806116                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2b                                                                                       |
| Conjugate           | DyLight 405                                                                                 |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                                                                                                                           |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                                                                     |
| Gene ID                 | 3066                                                                                                                                                                                                      |
| Gene Symbol             | HDAC2                                                                                                                                                                                                     |
| Species                 | Human, Rat                                                                                                                                                                                                |
| Specificity/Sensitivity | Detects human Histone Deacetylase 2/HDAC-2 in direct ELISAs and human and rat Histone Deacetylase 2/HDAC-2 Western blots. In direct ELISAs, no cross-reactivity with recombinant human HDAC1 is observed. |
| Immunogen               | <i>E. coli</i> -derived recombinant human Histone Deacetylase 2/HDAC-2 Pro386-Pro488<br>Accession # Q92769                                                                                                |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                         |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB7679E**

|                 |                                                           |
|-----------------|-----------------------------------------------------------|
| NBP2-27231V     | Mouse IgG2b Isotype Control (MPC-11) [DyLight 405]        |
| NBP1-72399-50ug | Recombinant Human Histone Deacetylase 2/HDAC2 His Protein |
| NBL1-11567      | H4/j Overexpression Lysate                                |
| NBL1-11480      | Histone Deacetylase 2/HDAC2 Overexpression Lysate         |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB7679E](http://www.novusbio.com/reviews/submit/FAB7679E)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

