

# Product Datasheet

## Protein Kinase A Regulatory Subunit II alpha Antibody (2394D) [DyLight 650] FAB8000W

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB8000W](http://www.novusbio.com/FAB8000W)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB8000W](http://www.novusbio.com/reviews/destination/FAB8000W)



**FAB8000W**

Protein Kinase A Regulatory Subunit II alpha Antibody (2394D) [DyLight 650]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	2394D
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 650
Purity	Protein A or G purified from cell culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Rabbit
Gene ID	5576
Species	Human, Mouse
Specificity/Sensitivity	Detects human PKA RII $\alpha$ in direct ELISAs. Detects human and mouse PKA RII $\alpha$ in Western blots.
Immunogen	Synthetic peptide containing Human/Mouse PKA RII $\alpha$ Accession # P13861
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunocytochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB8000W**

---

NBP2-24891C	Rabbit IgG Isotype Control [DyLight 650]
DYC1825-2	ERK1 [p Thr202, p Tyr204] [Biotin]
NBP1-48308	PKM2 Antibody - BSA Free
NBP1-97417	Vanilloid R1/TRPV1 Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB8000W](http://www.novusbio.com/reviews/submit/FAB8000W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

