

# Product Datasheet

## **B7-H7/HLA2 Antibody (907816) [FITC] FAB8084F**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB8084F](http://www.novusbio.com/FAB8084F)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB8084F](http://www.novusbio.com/reviews/destination/FAB8084F)



**FAB8084F**

B7-H7/HLA2 Antibody (907816) [FITC]

| Product Information         |                                                                                                                   |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                       |
| Storage                     | Store at 4C in the dark.                                                                                          |
| Clonality                   | Monoclonal                                                                                                        |
| Clone                       | 907816                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                |
| Isotype                     | IgG1                                                                                                              |
| Conjugate                   | FITC                                                                                                              |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                        |
| Buffer                      | PBS                                                                                                               |
| Product Description         |                                                                                                                   |
| Host                        | Mouse                                                                                                             |
| Gene ID                     | 11148                                                                                                             |
| Gene Symbol                 | HLA2                                                                                                              |
| Species                     | Human                                                                                                             |
| Specificity/Sensitivity     | Detects human B7-H7/HLA2 in direct ELISAs.                                                                        |
| Immunogen                   | Human embryonic kidney cell line HEK293-derived recombinant human B7-H7/HLA2<br>Met1-Asn344<br>Accession # Q9UM44 |
| Product Application Details |                                                                                                                   |
| Applications                | Immunocytochemistry                                                                                               |
| Recommended Dilutions       | Immunocytochemistry                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB8084F**

---

|               |                                         |
|---------------|-----------------------------------------|
| NBP1-96849    | Mouse IgG1 Isotype Control (MG1) [FITC] |
| NBP2-49187PEP | B7-H7/HLA2 Recombinant Protein Antigen  |
| 1849-NK-025   | NKp30/NCR3 [Unconjugated]               |
| 10109-B7-050  | B7-H7/HLA2 [Unconjugated]               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB8084F](http://www.novusbio.com/reviews/submit/FAB8084F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

