

# Product Datasheet

## **CXCR2/IL-8RB Antibody (866614R) [DyLight 650] FAB8110RW**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB8110RW](http://www.novusbio.com/FAB8110RW)

Updated 10/7/2024 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB8110RW](http://www.novusbio.com/reviews/destination/FAB8110RW)



**FAB8110RW**

CXCR2/IL-8RB Antibody (866614R) [DyLight 650]

| Product Information         |                                                                                                                                                                  |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                      |
| Storage                     | Store at 4C in the dark.                                                                                                                                         |
| Clonality                   | Monoclonal                                                                                                                                                       |
| Clone                       | 866614R                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                               |
| Isotype                     | IgG2b                                                                                                                                                            |
| Conjugate                   | DyLight 650                                                                                                                                                      |
| Purity                      | Protein A or G purified from cell culture supernatant                                                                                                            |
| Buffer                      | 50mM Sodium Borate                                                                                                                                               |
| Product Description         |                                                                                                                                                                  |
| Host                        | Mouse                                                                                                                                                            |
| Gene ID                     | 3579                                                                                                                                                             |
| Species                     | Rat                                                                                                                                                              |
| Specificity/Sensitivity     | Detects rat CXCR2/IL-8 RB in direct ELISAs. In Flow Cytometry, stains HEK293 cells transfected with rat CXCR2/IL-8 RB, but does not stain non-transfected cells. |
| Immunogen                   | NS0 mouse myeloma cell line transfected with rat CXCR2/IL-8 RB<br>Accession # P35407                                                                             |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                |
| Product Application Details |                                                                                                                                                                  |
| Applications                | Flow Cytometry, CyTOF-ready                                                                                                                                      |
| Recommended Dilutions       | Flow Cytometry, CyTOF-ready                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                           |





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to FAB8110RW

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| NBP2-27231C        | Mouse IgG2b Isotype Control (MPC-11) [DyLight 650]  |
| H00003579-P01-10ug | Recombinant Human CXCR2/IL-8RB GST (N-Term) Protein |
| 210-TA-005         | TNF-alpha [Unconjugated]                            |
| NBP3-06771         | Mouse CXCR2/IL-8RB ELISA Kit (Colorimetric)         |

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB8110RW](http://www.novusbio.com/reviews/submit/FAB8110RW)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

