

# Product Datasheet

## Flt-3/FIk-2/CD135 Antibody (66903) [DyLight 350] FAB812D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB812D](http://www.novusbio.com/FAB812D)

Updated 8/2/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB812D](http://www.novusbio.com/reviews/destination/FAB812D)



**FAB812D**

Flt-3/Flk-2/CD135 Antibody (66903) [DyLight 350]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	66903
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 350
Purity	Protein A or G purified from ascites
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	2322
Gene Symbol	FLT3
Species	Human
Specificity/Sensitivity	Detect human Flt $\alpha$ 3/Flk $\alpha$ 2 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant mouse Flt-3, recombinant human (rh) SCF R, rhPDGF R $\alpha$ , or rhPDGF R $\beta$ is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant human Flt $\alpha$ 3/Flk $\alpha$ 2 Asn27-Asn541 Accession # AAA18947
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Flow Cytometry, CyTOF-ready
Recommended Dilutions	Flow Cytometry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB812D**

---

NBP1-97005UV	Mouse IgG1 Isotype Control (MG1) [DyLight 350]
NBP2-57187PEP	FIt-3/Flk-2/CD135 Recombinant Protein Antigen
203-IL-010	IL-3 [Unconjugated]
368-ST-050/CF	FIt-3/Flk-2/CD135 [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB812D](http://www.novusbio.com/reviews/submit/FAB812D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

