

# Product Datasheet

## IL-10 Antibody (2472A) [DyLight 680] FAB91842Y

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB91842Y](http://www.novusbio.com/FAB91842Y)

Updated 5/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB91842Y](http://www.novusbio.com/reviews/destination/FAB91842Y)



**FAB91842Y**

IL-10 Antibody (2472A) [DyLight 680]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 2472A                                                                                       |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Conjugate           | DyLight 680                                                                                 |
| Purity              | Protein A or G purified from cell culture supernatant                                       |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                      |
|-------------------------|------------------------------------------------------------------------------------------------------|
| Host                    | Rabbit                                                                                               |
| Gene ID                 | 3586                                                                                                 |
| Species                 | Human                                                                                                |
| Specificity/Sensitivity | Detects human IL-10 in direct ELISAs.                                                                |
| Immunogen               | Human embryonic kidney cell, HEK293-derived human IL-10 protein<br>Met1-Asn178<br>Accession # P22301 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                    |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry                                     |
| Recommended Dilutions       | Western Blot, Immunohistochemistry                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB91842Y**

---

|                 |                                          |
|-----------------|------------------------------------------|
| NBP2-24891FR    | Rabbit IgG Isotype Control [DyLight 680] |
| NBP2-35101-10ug | Recombinant Mouse IL-10 Protein          |
| 210-TA-005      | TNF-alpha [Unconjugated]                 |
| DY417-05        | IL-10 [Biotin]                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB91842Y](http://www.novusbio.com/reviews/submit/FAB91842Y)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

