

Product Datasheet

ARG2 Antibody (3D10) - Azide and BSA Free H00000384-M05

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00000384-M05

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00000384-M05



H00000384-M05


ARG2 Antibody (3D10) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3D10
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	384
Gene Symbol	ARG2
Species	Human
Specificity/Sensitivity	ARG2 - arginase, type II (3D10)
Immunogen	ARG2 (AAH29050, 1 a.a. ~ 354 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MSLRGSLSRLLQTRVHSILKKS VHSVAVIGAPFSQQQKRKGVEHGPA AireAGL MKRLSSLGCHLKDFGDL SFTPVPKDDL YNNLIVNPRSVGLANQELAEVVSRAV SDGYSCVTLGGDHSLAIGTISGHARHCPDLCVWVDAHADINTPLTTSSGNLH GQPVSFLLRELQDKVPQLPGFSWIKPCISSASIVYIGLRDVPPEHFILKNYDIQY FSMRDIDRLGIQKVMERTFDLLIGKRQRPIHLSFDIDAFDPTLAPATGTPVVGGL TYREGMYIAEEIHNTGLLSALDLVEVNPQLATSEEEAKTTANLAVDVIASSFGQT REGGHIVYDQLPTPSSPDESENQARVRI
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	ELISA, Immunoprecipitation, Sandwich ELISA
Recommended Dilutions	ELISA, Immunoprecipitation, Sandwich ELISA
Application Notes	It has been used for ELISA.




Images

Immunoprecipitation: ARG2 Antibody (3D10) [H00000384-M05] - Analysis of ARG2 transfected lysate using anti-ARG2 monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with ARG2 rabbit polyclonal antibody.



Sandwich ELISA: ARG2 Antibody (3D10) [H00000384-M05] - Detection limit for recombinant GST tagged ARG2 is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00000384-M05

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP1-45256-0.25mg	Recombinant Human ARG2 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00000384-M05

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

