

Product Datasheet

DAZ1 Antibody (3A4) - Azide and BSA Free H00001617-M09

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001617-M09

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001617-M09




H00001617-M09

DAZ1 Antibody (3A4) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3A4
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	1617
Gene Symbol	DAZ1
Species	Human
Specificity/Sensitivity	DAZ1 - deleted in azoospermia 1
Immunogen	DAZ1 (AAH18119, 21 a.a. ~ 120 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. SSSAAASQGWVLPEGKIVPNTVFVGGIDARMDTEIGSCFGRYGSVKEVKIITN RTGVSKGYGFVSFVNDVDVQKIVGSQIHFGKKLKLGPAIRKQKLC
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	Antibody reactivity against tissue lysate and recombinant protein for WB. It has been used for IHC-P and ELISA.

Images


Western Blot: DAZ1 Antibody (3A4) [H00001617-M09] - DAZ1 monoclonal antibody (M09), clone 3A4. Western Blot analysis of DAZ1 expression in human placenta.



Immunohistochemistry-Paraffin: DAZ1 Antibody (3A4) [H00001617-M09] - Analysis of monoclonal antibody to DAZ1 on formalin-fixed paraffin-embedded human testis. Antibody concentration 3 ug/ml



ELISA: DAZ1 Antibody (3A4) [H00001617-M09] - Detection limit for recombinant GST tagged DAZ1 is approximately 1ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Error: Subreport could not be shown.

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001617-M09

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

