

Product Datasheet

Recombinant Human GABA-B R1 Protein H00002550-G01

Unit Size: 10 ug

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00002550-G01

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00002550-G01



H00002550-G01

Recombinant Human GABA-B R1 Protein

Product Information	
Unit Size	10 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -80C. Avoid freeze-thaw cycles.
Preservative	Glycerol
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	25 mM Tris-HCl pH8.0 in 2% glycerol.
Target Molecular Weight	95.1 kDa

Product Description	
Description	<p>An untagged recombinant protein corresponding to the amino acid sequence of (NP_068703.1) for Human GABA-B R1</p> <p>Source: <i>Wheat Germ (in vitro) with proprietary liposome technology</i></p> <p>Amino Acid Sequence: MGPGAPFARVGVWPLPLLVVMAAGVAPVWASHSPHLPRPHSRVPPHPSSERR AVYIGALFPMMSGGWPGGQACQPAVEMALEDVNSRRDILPDYELKLIHHSKCD PGQATKYLYELLYNDPIKIILMPGCSSVSTLVAEAARMWNLIVLSYGSSSPALSN RQRFPTFFRTHPSATLHNPTRVKLFEKWGWKKIATIQQTTEVFTSTLDDLEERV KEAGIEITFRQSFFSDPAVPVKNLKRQDARIIVGLFYETEARKVFCEVYKERLFG KKYVWFLIGWYADNWFKIYDPSINCTVDEMTEAVEGHITTEIVMLNPANTRSISN MTSQEFVEKLTKRLKRHPEETGGFQEAPLAYDAIWALALALNKTSGGGGRSGV RLEDFNYNNQTITDQIYRAMNSSFEGVSGHVVFDAASGRMAWTLIEQLQGGG YKKIGYYDSTKDDLSSKTDKWIGGSPPADQTLVIKTRFRLSQKLFISVSVLSSL GIVLAVVCLSFNIYNHVRVYIQNSQPNLNNLTAVGCSLALAAVFPLGLDGYHIGR NQFPFVCQARLWLLGLGFSLGYSMFTKIWWVHTVFTKKEEKKEWRKTLEPW KLYATVGLLVGMDVLTLAIWQIVDPLHRTIETFAKEEPKEDIDVSILPQLEHCSSR KMNTWLGIFYGYKGLLLLLGIFLAYETKSVSTEKINDHRAVGMAIYNVAVLCLITA PVTMILSSQQDAAFASLAIVFSSYITLVVLFVPMRRLITRGEWQSEAQDTMK TGSSTNNNEEEKSRLEKENRELEKIIAEKEERVSELRHQLQSRQQLRSRRHPP TPPEPSGGLPRGPPEPPDRLSCDGSRVHLLYK</p>
Gene ID	2550
Gene Symbol	GABBR1
Species	Human
Preparation Method	in vitro wheat germ expression system with proprietary liposome technology
Details of Functionality	This protein was produced in an in vitro wheat germ expression system that should preserve correct conformational folding that is necessary for biological function. While it is possible that this protein could display some level of activity, the functionality of this protein has not been explicitly measured or validated.
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Immunoaffinity Purification
Recommended Dilutions	Immunoaffinity Purification





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00002550-G01

H00002550-Q01-10ug	Recombinant Human GABA-B R1 GST (N-Term) Protein
203-IL-010	IL-3 [Unconjugated]
H00002550-M01	GABA-B R1 Antibody (2D7)
NB300-605	iNOS Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00002550-G01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

