

Product Datasheet

MOG 293T Cell Transient Overexpression Lysate H00004340-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00004340-T01

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00004340-T01



H00004340-T01**MOG 293T Cell Transient Overexpression Lysate****Product Information**

Unit Size	0.1 ml
Concentration	mg/ml
Storage	Store at -80C. Avoid freeze-thaw cycles.
Buffer	1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue)

Product Description

Description	Quality control test: Transient overexpression cell lysate was tested with Anti-MOG antibody by Western Blots. Plasmid: pCMV-MOG full-length
Gene ID	4340
Gene Symbol	MOG
Species	Human
Specificity/Sensitivity	MOG 293T Cell Transient Overexpression Lysate(Denatured)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Lysate Type	293T Cell Transient Overexpression
Lysate Protein State	Denatured

Product Application Details


Applications	Western Blot
Recommended Dilutions	Western Blot

Images

MOG transfected lysate



Lane 1: MOG transfected lysate (32.56 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00004340-T01

NBP2-46634PEP	MOG Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
H00004340-B01P	MOG Antibody
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00004340-T01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

