

# Product Datasheet

## Myosin Phosphatase 2 Antibody - Azide and BSA Free H00004660-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00004660-D01P](http://www.novusbio.com/H00004660-D01P)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00004660-D01P](http://www.novusbio.com/reviews/destination/H00004660-D01P)



**H00004660-D01P**

Myosin Phosphatase 2 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Rabbit
Gene ID	4660
Gene Symbol	PPP1R12B
Species	Human
Immunogen	PPP1R12B (ENSP00000349206, 1 a.a. - 386 a.a.) full-length human protein. MAELEHLGGKRAESARMRRAEQLRRWRGSLTEQEPAERRGAGRQPLTRRGS PRVRFEDGAVFLAACSSGDTDEVKLLARGADINTVNVVDGLTALHQACIDENLD MVKFLVENRANVNQQDNEGWTPLHAAASCGYLNIAEYFINHGASVGVNSEGE VPSDLAEEPAMKDLLLEQVKKQGVLDLEQSRKEEQMLQDARQWLNSGKIED VRQARSGATALHVAAAKGYSEVLRLLIQAGYELNVQDYDGWTPHAAAHWGV KEACSLAEALCDMDIRNKLQTPFDVADEGLVEHLELLQKKQNVLRSEKETRN KLIESDLNSKIQSGFFKNKEKMLYEEETPKSQEMEEENKESSSSSSEEEEGEDE ASESETEKEAVLFWPF
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	It has been used for WB and Functional.




**Images**

Western Blot: Myosin Phosphatase 2 Antibody [H00004660-D01P] - Analysis of PPP1R12B expression in human kidney.



Western Blot: Myosin Phosphatase 2 Antibody [H00004660-D01P] - Analysis of PPP1R12B expression in transfected 293T cell line by PPP1R12B polyclonal antibody. Lane 1: PPP1R12B transfected lysate (43.30 KDa). Lane 2: Non-transfected lysate.





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to H00004660-D01P

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00004660-P01-10ug	Recombinant Human Myosin Phosphatase 2 GST (N-Term) Protein

---

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00004660-D01P](http://www.novusbio.com/reviews/submit/H00004660-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

