

Product Datasheet

PLRP1 Antibody - Azide and BSA Free H00005407-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00005407-B01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00005407-B01P




H00005407-B01P

PLRP1 Antibody - Azide and BSA Free


Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Product Description	
Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Mouse
Gene ID	5407
Gene Symbol	PNLIPRP1
Species	Human
Specificity/Sensitivity	PNLIPRP1 - pancreatic lipase-related protein 1,
Immunogen	PNLIPRP1 (AAH25784.1, 1 a.a. - 467 a.a.) full-length human protein. MLIFWTITLFLLGAAKGKEVCYEDLGCFSDTEPWGGTAIRPLKILPWSPEKIGTR FLLYTNENPNNFQILLSDPSTIEASNFQMDRKRFIHGFIDKGDSEWVTDMCK KLFEVEEVNVCICVDWKKGSQATYTQAANNVRVVGQAQVAQMLDILLTEYSYPPS KVHLIGHSLGAHVAGEAGSKTPGLSRITGLDPVEASFESTPEEVRLDPSDADFV DVIHTDAAPLIPFLGFGTNQQMGHLDFPNGGESMPGCKKNALSQIVDLGIW AGTRDFVACNHLRSYKYLESILNPDGFAAYPCTSYKSFESDKCFPCPDQGCP QMGHYADKFAGRTSEEQQKFFLNTGEASN FARWRYGVSITLSGRTATGQIKVA LFGNKGNTHQYSIFRGILKPGSTHSYEFDAKLDVGTIDKVKFLWNNNVINPTLPK VGATKITVQKGEEKTVYNFCSEDTVREDTLLTLTPC
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Antibody reactivity against transfected lysate in WB. It has been used for ELISA.

Images

Western Blot: PLRP1 Antibody [H00005407-B01P] - Analysis of PNLIPRP1 expression in human pancreas.



Western Blot: PLRP1 Antibody [H00005407-B01P] - Analysis of PNLIPRP1 expression in transfected 293T cell line by PNLIPRP1 polyclonal antibody. Lane 1: PNLIPRP1 transfected lysate(51.37 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00005407-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBL1-14544	PLRP1 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00005407-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

