

Product Datasheet

PTMS Antibody (3H5) - Azide and BSA Free H00005763-M17-100ug

Unit Size: 100 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00005763-M17-100ug

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00005763-M17-100ug



H00005763-M17-100ug

PTMS Antibody (3H5) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3H5
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	In 1x PBS, pH 7.4


Product Description	
Description	Novus Biologicals Mouse PTMS Antibody (3H5) - Azide and BSA Free (H00005763-M17-100ug) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5763
Gene Symbol	PTMS
Species	Human
Immunogen	PTMS (NP_002815.3, 1 a.a. ~ 102 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MSEKSVEAAAELSAKDLKEKKEKVEEKASRKERKKEVVEEEEENGAEIIIIIIETA EDGEEDDEGEEEDDEEEEDDEDEGPALKRAAEEEDDEADPKRQKTENGASA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence




Images


Western Blot: PTMS Antibody (3H5) [H00005763-M17-100ug] - Detection against Immunogen (37.9 KDa) .




Immunocytochemistry/Immunofluorescence: PTMS Antibody (3H5) [H00005763-M17-100ug] - Monoclonal antibody to PTMS on HeLa cell . [antibody concentration 10 ug/ml]



Western Blot: PTMS Antibody (3H5) [H00005763-M17-100ug] - Analysis of PTMS expression in transfected 293T cell line by PTMS monoclonal antibody (M17), clone 3H5. Lane 1: PTMS transfected lysate (Predicted MW: 11.5 KDa). Lane 2: Non-transfected lysate.



ELISA: PTMS Antibody (3H5) [H00005763-M17-100ug] - Detection limit for recombinant GST tagged PTMS is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00005763-M17-100ug

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00005763-M17-100ug

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

