

# Product Datasheet

## RARRES3 Antibody (1H5) - Azide and BSA Free H00005920-M10

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00005920-M10](http://www.novusbio.com/H00005920-M10)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00005920-M10](http://www.novusbio.com/reviews/destination/H00005920-M10)



**H00005920-M10**

RARRES3 Antibody (1H5) - Azide and BSA Free


| Product Information |  |
|---------------------|--|
| Unit Size           | 0.1 mg   |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.   |
| Clonality           | Monoclonal   |
| Clone               | 1H5  |
| Preservative        | No Preservative  |
| Isotype             | IgG2a Kappa  |
| Purity              | IgG purified   |
| Buffer              | In 1x PBS, pH 7.4  |

| Product Description     |   |
|-------------------------|---|
| Description             | Quality control test: Antibody Reactive Against Recombinant Protein.  |
| Host                    | Mouse   |
| Gene ID                 | 5920  |
| Gene Symbol             | RARRES3   |
| Species                 | Human   |
| Specificity/Sensitivity | RARRES3 - retinoic acid receptor responder (tazarotene induced) 3 (1H5)   |
| Immunogen               | RARRES3 (AAH09678, 1 a.a. ~ 164 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>MASPHQEPKPGDLIEIFRLGYEHWALYIGDGYVIHLAPPSEYPGAGSSSVFSVL<br>SNSAEVKRERLEDVVGCCYRVNNSLDHEYQPRPVEVISSAKEMVGQKMKY<br>SIVSRNCEHFVTQLRYGKSRCKQVEKAKVEVGVATALGILVVAGCSFAIRRYQK<br>KATA |
| Notes                   | This product is produced by and distributed for Abnova, a company based in Taiwan.  |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, ELISA  |
| Recommended Dilutions       | Western Blot 1:500, ELISA  |
| Application Notes           | Antibody reactivity against Recombinant Protein with GST tag on ELISA and WB. GST tag alone is used as a negative control. |

**Images**

Western Blot: RARRES3 Antibody (1H5) [H00005920-M10] - Analysis of RARRES3 expression in transfected 293T cell line by RARRES3 monoclonal antibody (M10), clone 1H5. Lane 1: RARRES3 transfected lysate (18.2 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

[nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

[canada.inquires@bio-techne.com](mailto:canada.inquires@bio-techne.com)

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

[info.EMEA@bio-techne.com](mailto:info.EMEA@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)

Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)

Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)

General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to H00005920-M10**

---

|                  |  |
|------------------|--|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]       |
| NB7539           | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK)           |
| 210-TA-005       | TNF-alpha [Unconjugated]                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00005920-M10](http://www.novusbio.com/reviews/submit/H00005920-M10)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

