

Product Datasheet

RHD Antibody - Azide and BSA Free H00006007-B01P-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00006007-B01P-50ug

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00006007-B01P-50ug



H00006007-B01P-50ug

RHD Antibody - Azide and BSA Free

Product Information

Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.4)

Product Description


Description	Novus Biologicals Mouse RHD Antibody - Azide and BSA Free (H00006007-B01P-50ug) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	6007
Gene Symbol	RHD
Species	Human
Immunogen	RHD (AAI60075.1, 1 a.a. - 417 a.a.) full-length human protein. MSSKYPRSVRRCLPLWALTLEAAILLFYFFTHYDASLEDQKGLVASVYQVGQDL TVMAAIGLGFLTSSFRRHSWSSVAFNLFMLALGVQWAILLDGFLSQFPSGKVV TLFSIRLATMSALSVLISVDAVLGKVNLAQLVVMVLVEVTALGNLRMVISNIFNTD YHMNMMHIYVFAAYFGLSVAWCLPKPLPEGTEKDQQTATIPSLSAMLGALFLW MFWPSFNSALLRSPIERKNAVFNTYYAVAVSVVTAISGSSLAHPQGKISKTYVH SAVLAGGVAVGTSCHLIPSPWLAMVLGLVAGLISVGGAKYLPGCCNRVLGIPHS SIMGYNFSLGLLGEIYIVLLVLDTVGAGNGMIGFQVLLSIGELSLAIVIALTSGLL TGLLLNLKIWKAPHEAKYFDDQVFWKFPHLAVGF
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot

Images

Western Blot: RHD Antibody [H00006007-B01P-50ug] - Analysis of RHD expression in transfected 293T cell line (H00006007-T02) by RHD MaxPab polyclonal antibody. Lane 1: RHD transfected lysate(45.87 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00006007-B01P-50ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00006007-B01P-50ug

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



