

Product Datasheet

SKIV2L Antibody (1E5) - Azide and BSA Free H00006499-M05-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00006499-M05-50ug

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00006499-M05-50ug



H00006499-M05-50ug


SKIV2L Antibody (1E5) - Azide and BSA Free

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1E5
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	6499
Gene Symbol	SKIV2L
Species	Human
Immunogen	SKIV2L (NP_008860, 1125 a.a. ~ 1233 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. DQLPNTLKQGIERVRAVAKRIGEVQVACGLNQTVEEFVVGELNFGFLVEVVYEWARGMPSSELAGLSGTPEGLVVRCIQRLAEMCRSLRGAARLVGEPVLGAKMETA ATLL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot, ELISA




Images


Western Blot: SKIV2L Antibody (1E5) [H00006499-M05-50ug] - Detection against Immunogen (37.73 KDa) .




Western Blot: SKIV2L Antibody (1E5) [H00006499-M05-50ug] - Analysis of SKIV2L expression in transfected 293T cell line by SKIV2L monoclonal antibody (M05), clone 1E5. Lane 1: SKIV2L transfected lysate (Predicted MW: 137.06 KDa). Lane 2: Non-transfected lysate.



Western Blot: SKIV2L Antibody (1E5) [H00006499-M05-50ug] - SKIV2L monoclonal antibody (M05), clone 1E5. Western Blot analysis of SKIV2L expression in Hela S3 NE.



ELISA: SKIV2L Antibody (1E5) [H00006499-M05-50ug] - Detection limit for recombinant GST tagged SKIV2L is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00006499-M05-50ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
H00006499-Q01-10ug	Recombinant Human SKIV2L GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00006499-M05-50ug

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

