

Product Datasheet

STAT4 Antibody (S1) H00006775-M02A-200ul

Unit Size: 200 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00006775-M02A-200ul

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00006775-M02A-200ul



H00006775-M02A-200ul**STAT4 Antibody (S1)**

Product Information	
Unit Size	200 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S1
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In ascites fluid


Product Description	
Host	Mouse
Gene ID	6775
Gene Symbol	STAT4
Species	Human
Specificity/Sensitivity	STAT4 - signal transducer and activator of transcription 4 (2H4-1G8)
Immunogen	STAT4 (AAH31212.1, 1 a.a. ~ 748 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MSQWNQVQQLEIKFLEQVDQFYDDNFPMEIRHLLAQWIENQDWEAASNNETM ATILLQNLLIQLDEQLGRVSKEKNLLLIHNLKRIRKVLQGKFHGNPMHVAVVISNC LREERRILAAANMPVQGPLEKSLQSSSVSERQRNVEHKVAAIKNSVQMTEQDT KYLEDLQDEFDYYRYKTIQTMDQSDKNSAMVNQEVLTLQEMLSLDFKRKEALS KMTQIIHETDLLMNTMLIEELQDWKRRQQIACIGGPLHNGLDQLQNCFTLLAESL FQLRRQLEKLEEQSTKMTYEGDPIPMQRTHMLERVTFIYNLFKNSFVVERQPC MPTHPQRPLVLKTLIQFTVKLRLLIKLPELNYQVKVKASIDKNVSTLSNRRFVLC GTNVKAMSIEESSNGSLSVFRHLQPKEMKSSAGGKGNEGCHMVTEELHSITF ETQICLYGLTIDLETSSLPVVMISNVSQLPNAWASIIWYNVSTNDSQNLVFFNNP PPATLSQLLEVMSWQFSSYVGRGLNSDQLHMLAEKLTQSSYSDGHILTWAKF CKEHLPGKSFTFTWLEAILDLIKKHILPLWIDGYVMGFVSKEKERLLLKDKMPG TFLLRFSESHLGGITFTWVDHSESGEVRFHVSVEPYNKGRLSALPFADILRDYKVI MAENIPENPLKYLYPDIPKDKAFGKHYSQPCVSRPTERGDKGYPVSVFIPIST IRSDSTEPHSPDLLPMSPSVYAVLRENLSPTTIETAMKSPYSAE
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot, ELISA



Images

Western Blot: STAT4 Antibody (S1) [H00006775-M02A-200ul] -
Detection against Immunogen (108.02 KDa) .





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00006775-M02A-200ul

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00006775-P01-10ug	Recombinant Human STAT4 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00006775-M02A-200ul

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

