

Product Datasheet

SULT1A3 Antibody - Azide and BSA Free H00006818-B02P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00006818-B02P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00006818-B02P



H00006818-B02P

SULT1A3 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Mouse
Gene ID	6818
Gene Symbol	SULT1A3
Species	Human
Immunogen	SULT1A3 (NP_003157.1, 1 a.a. - 295 a.a.) full-length human protein. MELIQDTSRPPLEYVKGVPLIKYFAEALGPLQSFQARPDDLINTYPKSGTTWV SQILDMIYQGGDLEKCNRAPIYVRVPFLEVNDPGEPSGLETLKDTPPPRLIKSHL PLALLPQTLLDQKVKVYVARNPKDVAVSYYHFHRMEKAHPEPGTWDSFLEKF MAGEVSYGSWYQHVVQEWELSRTHPVLVLYFYEDMKENPKREIQKILEFVGRS LPEETMDFMVQHTSFKEMKKNPMTNYTTVPQELMDHSISPFMRKGMAGDWK TTFTVAQNERFDADYAEKMAGCSLSFRSEL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Antibody reactivity against transfected lysate for WB and recombinant protein with GST tag for ELISA and WB. GST tag alone is used as a negative control.



Images


Western Blot: SULT1A3 Antibody [H00006818-B02P] - Analysis of SULT1A3 expression in human liver.



Western Blot: SULT1A3 Antibody [H00006818-B02P] - Analysis of SULT1A3 expression in HepG2.



Western Blot: SULT1A3 Antibody [H00006818-B02P] - Analysis of SULT1A3 expression in transfected 293T cell line by SULT1A3 polyclonal antibody. Lane 1: SULT1A3 transfected lysate(32.45 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00006818-B02P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
7027-GT-020	ST8 alpha-2,8-Sialyltransferase 4/ST8SIA4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00006818-B02P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

