

Product Datasheet

ZSCAN21/ZFP38 Antibody (4B11) - Azide and BSA Free H00007589-M08

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00007589-M08

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00007589-M08



H00007589-M08


ZSCAN21/ZFP38 Antibody (4B11) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4B11
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	7589
Gene Symbol	ZSCAN21
Species	Human
Specificity/Sensitivity	ZNF38 (4B11)
Immunogen	ZSCAN21 (AAH47309.1, 1 a.a. ~ 473 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MMTKVLGMAPVLGPRPPQEQVGPLMVKVEEKEEKGYLPSLEMFRQRFQF GYHDTGPGPREALSQLRVLCCEWLRPEIHTKEQILELLVLEQFLTILPQELQAWV QEHCPESAEEAVTLLEDLERELDEPGHQVSTPPNEQKPVWEKISSSGTAKESP SSMQPQPLETSHKYESWGPLYIQESGEEQEFAQDPRKVRDCRLSTQHEESAD EQKGSEAEGLKGDIIISVIIANKPEASLERQCVNLENEKGTKPPLQEAGSKKGRE SVPTKPTPGERRYICAECGKAFSNSSNLTKHRRHTHTGEKPYVCTKCGKAFSHS SNLTLHYRTHLVDRPYDCKCGKAFGQSSDLLKHQRMHTEEAPYQCKDCGKAF SGKGLIRHYRIHTGEKPYQCNECGKSFSQHAGLSSHQRLHTGEKPYKCKEC GKA FNHSSNFNKHRIHTGEKPYWCHHCGKTFCSKSNLSKHQRVHTGEGEAP
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot, ELISA
Application Notes	Antibody reactivity against recombinant protein on ELISA.




Images

Western Blot: ZSCAN21/ZFP38 Antibody (4B11) [H00007589-M08] - Analysis of ZSCAN21 expression in transfected 293T cell line by ZNF38 monoclonal antibody (M08), clone 4B11. Lane 1: ZSCAN21 transfected lysate (Predicted MW: 53.7 KDa). Lane 2: Non-transfected lysate.



ELISA: ZSCAN21/ZFP38 Antibody (4B11) [H00007589-M08] - Detection limit for recombinant GST tagged ZSCAN21 is 0.1 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00007589-M08

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP1-84182PEP	ZSCAN21/ZFP38 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00007589-M08

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

