

Product Datasheet

SNAP23 Antibody (2F5-3D4) - Azide and BSA Free H00008773-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00008773-M01

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00008773-M01



H00008773-M01


SNAP23 Antibody (2F5-3D4) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2F5-3D4
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	8773
Gene Symbol	SNAP23
Species	Human
Specificity/Sensitivity	SNAP23 - synaptosomal-associated protein, 23kDa
Immunogen	SNAP23 (AAH00148, 1 a.a. ~ 211 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MDNLSSEIIQQR AHQITDESLESTRILGLAIESQDAGTKTITMLDEQKEQLNRI EEGLDQINKDMRETEKLTLELNKCCGLCVCPCNRTKNFESGKAYKTTWGDGG ENSPCNVVSKQPGPVTNGQLQQPTTGAASGGYIKRITNDAREDEMEENLTQV GSILGNLKD MALNIGNEIDAQNPQIKRITDKADTNRDRIDIANARAKKLIDS
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, Sandwich ELISA
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Paraffin, Sandwich ELISA
Application Notes	Antibody reactivity against cell lysate and recombinant protein for WB. It has also been used for IHC-P and ELISA.




Images


Western Blot: SNAP23 Antibody (2F5-3D4) [H00008773-M01] - Analysis of SNAP23 expression in transfected 293T cell line by SNAP23 monoclonal antibody (M01), clone 2F5-3D4. Lane 1: SNAP23 transfected lysate (Predicted MW: 23.4 KDa). Lane 2: Non-transfected lysate.




Immunohistochemistry-Paraffin: SNAP23 Antibody (2F5-3D4) [H00008773-M01] - Analysis of monoclonal antibody to SNAP23 on formalin-fixed paraffin-embedded human spleen. Antibody concentration 3 ug/ml.




Western Blot: SNAP23 Antibody (2F5-3D4) [H00008773-M01] - SNAP23 expression in Raw 264.7.



Immunoprecipitation: SNAP23 Antibody (2F5-3D4) [H00008773-M01] - Analysis of SNAP23 transfected lysate using anti-SNAP23 monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with SNAP23 MaxPab rabbit polyclonal antibody.



Sandwich ELISA: SNAP23 Antibody (2F5-3D4) [H00008773-M01] - Detection limit for recombinant GST tagged SNAP23 is approximately 0.03ng/ml as a capture antibody.



Publications

Larkin C Examining the TRPV3 cation channel in a culture model of human epidermis Thesis 2023-01-01 (ICC/IF, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00008773-M01

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBC1-18347-0.05mg	Recombinant Human SNAP23 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00008773-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

