

# Product Datasheet

## GKAP/DLGAP1 Antibody (2A6) - Azide and BSA Free H00009229-M02-100ug

Unit Size: 100 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00009229-M02-100ug](http://www.novusbio.com/H00009229-M02-100ug)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00009229-M02-100ug](http://www.novusbio.com/reviews/destination/H00009229-M02-100ug)



**H00009229-M02-100ug**


GKAP/DLGAP1 Antibody (2A6) - Azide and BSA Free

| Product Information         |                                                                                                                                                                                                                                |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ug                                                                                                                                                                                                                         |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                                               |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                                                   |
| Clonality                   | Monoclonal                                                                                                                                                                                                                     |
| Clone                       | 2A6                                                                                                                                                                                                                            |
| Preservative                | No Preservative                                                                                                                                                                                                                |
| Isotype                     | IgG2a Kappa                                                                                                                                                                                                                    |
| Purity                      | Protein A or G purified                                                                                                                                                                                                        |
| Buffer                      | In 1x PBS, pH 7.4                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                |
| Host                        | Mouse                                                                                                                                                                                                                          |
| Gene ID                     | 9229                                                                                                                                                                                                                           |
| Gene Symbol                 | DLGAP1                                                                                                                                                                                                                         |
| Species                     | Human                                                                                                                                                                                                                          |
| Immunogen                   | DLGAP1 (NP_004737.2, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>MKGLSGSRSHHHGVTCDSACDSLHHSDRKPYLLSPVEHHPADHPYYTQRNS<br>FQAECVGPFSFDPLASSTFPRRHYSQQELKDECALVPRTLATKANRIP |
| Notes                       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                             |
| Product Application Details |                                                                                                                                                                                                                                |
| Applications                | Western Blot, ELISA                                                                                                                                                                                                            |
| Recommended Dilutions       | Western Blot, ELISA                                                                                                                                                                                                            |




## Images

Western Blot: GKAP/DLGAP1 Antibody (2A6) [H00009229-M02-100ug] - Detection against Immunogen (36.74 KDa) .



ELISA: GKAP/DLGAP1 Antibody (2A6) [H00009229-M02-100ug] - Detection limit for recombinant GST tagged DLGAP1 is 0.03 ng/ml as a capture antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00009229-M02-100ug**

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK)                |
| NBP1-76911PEP    | GKAP/DLGAP1 Antibody Blocking Peptide                   |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00009229-M02-100ug](http://www.novusbio.com/reviews/submit/H00009229-M02-100ug)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

