

Product Datasheet

Recombinant Human GPR37L1 Protein H00009283-G01

Unit Size: 10 ug

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009283-G01

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009283-G01



H00009283-G01

Recombinant Human GPR37L1 Protein

| Product Information | |
|-----------------------------|--|
| Unit Size | 10 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -80C. Avoid freeze-thaw cycles. |
| Preservative | Glycerol |
| Purity | >80% by SDS-PAGE and Coomassie blue staining |
| Buffer | 25 mM Tris-HCl of pH8.0 containing 2% glycerol. |
| Target Molecular Weight | 52.8 kDa |
| Product Description | |
| Description | <p>An untagged recombinant protein corresponding to the amino acid sequence of (AAH50334.1) for Human GPR37L1</p> <p>Source: <i>Wheat Germ (in vitro) with proprietary liposome technology</i></p> <p>Amino Acid Sequence: MRWLWPLAVSLAVILAVGLSRVSGGAPLHLGRHRAETQEQQSRSKRGTEDEE AKGVQQYVPEEWAHEYPRPIHPAGLQPTKPLVATSPNPGRDGGTPDSGQELRG NLTGAPGQRLQIQNPLYPVTESSYSAYAIMLLALVFAVGIVGNLSVMCIVWHS YYLKSAWNSILASLALWDFLVLFFCLPIVIFNEITKQRLLGDVSCRAVPFMEVSSL GVTTFSLCALGIDRFHVATSTLPKVRPIERCQSILAKLAVIWWGSMTLAVPELLL WQLAQEPAPTMGTLDSCIMKPSASLPESLYSLVMTYQNARMWWWYFGCYFCLP ILFTVTCQLVTWRVRGPPGRKSECRASKHEQCESQLNSTVVGLTVVYAFCTLP ENVCNIVVAYLSTELTRQTLDLLGLINQFSTFFKGAITPVLLLCICRPLGQAFLDC CCCCCEECGGASEASAANGSDNKLKTEVSSSIYFHKPRESPLLPLGTPC</p> |
| Gene ID | 9283 |
| Gene Symbol | GPR37L1 |
| Species | Human |
| Preparation Method | in vitro wheat germ expression system with proprietary liposome technology |
| Details of Functionality | This protein was produced in an in vitro wheat germ expression system that should preserve correct conformational folding that is necessary for biological function. While it is possible that this protein could display some level of activity, the functionality of this protein has not been explicitly measured or validated. |
| Notes | This product is produced by and distributed for Abnova, a company based in Taiwan. |
| Product Application Details | |
| Applications | Immunoaffinity Purification |
| Recommended Dilutions | Immunoaffinity Purification |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00009283-G01

| | |
|--------------------|--|
| H00009283-Q01-10ug | Recombinant Human GPR37L1 GST (N-Term) Protein |
| 233-FB-025 | FGF basic/FGF2/bFGF [Unconjugated] |
| MAB4449 | GPR37L1 Antibody (444532) [Unconjugated] |
| 314-BP-010 | BMP-4 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009283-G01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

