

Product Datasheet

GTF2IRD1 Antibody (1D2) - Azide and BSA Free H00009569-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009569-M01

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009569-M01




H00009569-M01**GTF2IRD1 Antibody (1D2) - Azide and BSA Free**


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1D2
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	9569
Gene Symbol	GTF2IRD1
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	GTF2IRD1 - GTF2I repeat domain containing 1
Immunogen	GTF2IRD1 (NP_005676, 431 a.a. ~ 530 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. IFTGNKFTKDTTKLEPASPPEDTSAEVSRA TVLDLAGNARSDK GSMSEDCGPG TSGELGGLRPIKIEPEDLDIIQVTPDPSP TSEEMTDSMPGHLPS E D
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody reactive against cell lysate, transfected lysate and recombinant protein for western blot. It has also been used for ELISA.

Images

Western Blot: GTF2IRD1 Antibody (1D2) [H00009569-M01] - GTF2IRD1 monoclonal antibody (M01), clone 1D2 Analysis of GTF2IRD1 expression in Jurkat.



Western Blot: GTF2IRD1 Antibody (1D2) [H00009569-M01] - Analysis of GTF2IRD1 expression in transfected 293T cell line by GTF2IRD1 monoclonal antibody (M01), clone 1D2. Lane 1: GTF2IRD1 transfected lysate(106 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00009569-M01

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
H00009569-P01-10ug	Recombinant Human GTF2IRD1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009569-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

