

# Product Datasheet

## PUM3 Antibody H00009933-B01

Unit Size: 0.05 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00009933-B01](http://www.novusbio.com/H00009933-B01)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00009933-B01](http://www.novusbio.com/reviews/destination/H00009933-B01)



**H00009933-B01**

## PUM3 Antibody


Product Information	
Unit Size	0.05 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	No additive

Product Description	
Description	Novus Biologicals Mouse PUM3 Antibody (H00009933-B01) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9933
Gene Symbol	KIAA0020
Species	Human
Specificity/Sensitivity	KIAA0020 - KIAA0020,
Immunogen	KIAA0020 (AAH16137.2, 1 a.a. - 648 a.a.) full-length human protein. MEVKGKKQFTGKNTKTAQEKNRFHKNSDSGSSKTFPTRKVAKEGGPKVTSRN FEKSITKLGKKGVKQFKNKQQGDKSPKNKFQANKFNKRRKFQPDGRSDESA AKKPKWDDFKKKKKELKQSRQLSDKTNYDIVVRAKQMWEILRRKDCDKEKRV KLMSDLQKLIQGGIKTIAFAHDSTRVIQCYIQYGNEEQRKQAFEELRDDLVLSK AKYSRNIVKKFLMYGSKPQIAEIIIRSFKGHVRKMLRHAEASAIVEYAYNDKAILE QRNMLTEELYGNTFQLYKSADHPTLDKVLLELQPEKLELIMDEMKGILTQMAQKE AVIKHSLVHKVFLDFFTYAPPKLRSEMIEAIREAVVYLAHTHDGARVAMHCLWH GTPKDRKIVVKTMYVEKVANGQYSHLVLLAAFDCIDDTKLVKQIIISEIISLPSI VNDKYGRKVLLYLLSPRDPAAHTVREIIEVLQKGDGNAHSSKDETEVRRRELLESI SPALLSYLQEAHAQEVVLDKSACVLVSDILGSATGDVQPTMNAIASLAATGLHPG GKDGEHIAEHPAGHLVLKWLIEQDKMKENGREGCFAKTLVEHVGMKNLKS WASVNRGAILSSLLQSCDLEVANKVKAALKSLIPTLEKTKSTSKGIEILLEKLST
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/ Immunofluorescence
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.

## Images

Western Blot: KIAA0020 Antibody [H00009933-B01] - Analysis of KIAA0020 expression in transfected 293T cell line by KIAA0020 polyclonal antibody. Lane 1: KIAA0020 transfected lysate(71.28 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: KIAA0020 Antibody [H00009933-B01] - Analysis of purified antibody to KIAA0020 on HeLa cell. (antibody concentration 10 ug/ml)



Western Blot: KIAA0020 Antibody [H00009933-B01] - Western Blot analysis of KIAA0020 expression in transfected 293T cell line by KIAA0020 MaxPab polyclonal antibody (B01). Lane 1: KIAA0020 transfected lysate (73.60 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00009933-B01**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97019-5mg	Mouse IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00009933-B01](http://www.novusbio.com/reviews/submit/H00009933-B01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

