

Product Datasheet

RASGRP2 Antibody (3D8) - Azide and BSA Free H00010235-M09

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010235-M09

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010235-M09




H00010235-M09**RASGRP2 Antibody (3D8) - Azide and BSA Free**

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3D8
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse RASGRP2 Antibody (3D8) - Azide and BSA Free (H00010235-M09) is a monoclonal antibody validated for use in WB and ELISA. Anti-RASGRP2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10235
Gene Symbol	RASGRP2
Species	Human, Canine
Specificity/Sensitivity	RASGRP2 (3D8)
Immunogen	RASGRP2 (NP_005816, 65 a.a. ~ 164 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. GTLDLDKGCTVEELLRGCI AFDDSGKVRDPQLVRMFLMMHPWYIPSSQLAAK LLHIYQQSRKDNSNSLQVKTCHLVRYWISAFPAEFDLNPELAEQIKE
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody reactive against recombinant protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.




Images

Western Blot: RASGRP2 Antibody (3D8) [H00010235-M09] - Analysis of RASGRP2 expression in transfected 293T cell line by RASGRP2 monoclonal antibody (M09), clone 3D8. Lane 1: RASGRP2 transfected lysate (75.5 KDa). Lane 2: Non-transfected lysate.



ELISA: RASGRP2 Antibody (3D8) [H00010235-M09] - Detection limit for recombinant GST tagged RASGRP2 is approximately 0.3ng/ml as a capture antibody.



Publications

Boudreaux MK, Catalano JL, Klok M. Calcium-Diacylglycerol Guanine Nucleotide Exchange Factor I (CaDAG-GEFI) gene mutations associated with loss of function in canine platelets. Transl Res;150(2):81-92. 2007-08-01 [PMID: 17656327]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00010235-M09

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP1-86661PEP	RASGRP2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010235-M09

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

