

Product Datasheet

Tubulin Beta 2C Antibody (1G3) - Azide and BSA Free H00010383-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010383-M02

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010383-M02




H00010383-M02

Tubulin Beta 2C Antibody (1G3) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1G3
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	10383
Gene Symbol	TUBB4B
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	TUBB2C - tubulin, beta, 2
Immunogen	TUBB4B (AAH01911, 1 a.a. ~ 445 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MREIVHLQAGQCGNQIGAKFWEVISDEHGIDPTGTYHGSDLQLERINVYYNEA TGGKYVPRAVLVDLEPGTMDSVRSRSGPFGQIFRPDNFVFGQSGAGNNWAKGH YTEGAELVDSVLDVVRKEAESCDCQLQGFQLTHSLGGGTGSGMGTLISKIREE YPDRIMNTFSVVPSPKVS DTVVEPYNATLSVHQLVENTDETYCIDNEALYDICFR TLKLTTPTYGDLNHLVSATMSGVTTCLRFPGLNADLRKLVNMMVFPRLHFF MPGFAPLTSRGSQQYRALTVP ELTQQMFD AKNMMAACDPRHGRYLTVAAVFR GRMSMKEVDEQMLNVQNKNS SYFVEWIPNNVKTAVCDIPRGLKMSATFIGNS TAIQELFKRISEQFTAMFRRKAFLHWYTGEGMDEMEFTEAESNMNDLVSEYQQ YQDATAEEEEGEFEEEEEEEEVA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500, ELISA, Immunocytochemistry/ Immunofluorescence
Application Notes	Antibody reactive against cell lysate and recombinant protein for Western Blot. Has also been used for immunofluorescence and ELISA.

Images

Western Blot: Tubulin Beta 2C Antibody (1G3) [H00010383-M02] - TUBB4B monoclonal antibody (M02), clone 1G3 Analysis of TUBB4B expression in HeLa.



Immunocytochemistry/Immunofluorescence: Tubulin Beta 2C Antibody (1G3) [H00010383-M02] - Analysis of monoclonal antibody to TUBB4B on HeLa cell. Antibody concentration 30 ug/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00010383-M02

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00010383-P01-10ug	Recombinant Human Tubulin Beta 2C GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010383-M02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

