

Product Datasheet

SUGT1 Antibody - Azide and BSA Free H00010910-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010910-D01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010910-D01P



H00010910-D01P

SUGT1 Antibody - Azide and BSA Free

Product Information

Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description


Host	Rabbit
Gene ID	10910
Gene Symbol	SUGT1
Species	Human
Specificity/Sensitivity	SUGT1 - SGT1, suppressor of G2 allele of SKP1 (<i>S. cerevisiae</i>),
Immunogen	SUGT1 (NP_006695.1, 1 a.a. - 333 a.a.) full-length human protein. MAAAAAGTATSQRFFQSFS DALIDEDPQA ALEELTKALEQKPDDAQYYCQRAY CHILLGNYCVAVADAKKSLELNPNNSTAMLRKGICEYHEKNYAAALETFTGQK LDSADANFSVWIKRCQEAQNGSESEVWTHQSKIKYDWYQTESQVVITLMIKNV QKNDVNVEFSEKELSALVKLPSGEDYNLKLELLHPIPEQSTFKVLSTKIEIKLKK PEAVRWEKLEGQGQDVPTPKQFVADVKNLYPSSSPYTRNWDKLVGEIKEEEKN EKLEGDAALNRLFQQIYSDGSDEVKRAMNKSFMESGGTVLSTNWSVDVGKRKV EINPPDDMEWKKY
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1:500
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.

Images

Western Blot: SUGT1 Antibody [H00010910-D01P] - Analysis of SUGT1 expression in transfected 293T cell line by SUGT1 polyclonal antibody. Lane 1: SUGT1 transfected lysate (37.80 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to H00010910-D01P

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-50979-0.1mg	Recombinant Human SUGT1 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010910-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

